

## **SPT-1 Sensor Pad Transmitter**

## Giving You a Peace of Mind



- Compatible with bed, chair, and floor mat sensor systems
- Detects when patient gets up from the bed/chair pad, or steps on floor mat
- Connect type: handset jack (RJ-9)

The SPT-1 is a wireless sensor pad transmitter that can be connected to a bed/chair sensor pad or alarming floor mat to create a bed, chair, and wander detection wireless alarm system. When a resident or patient gets up from the bed/chair pad, or pressure is added to a floor mat, a wireless signal and alarm will go off and a caregiver would be immediately alerted of the situation. The SPT-1 can help track and monitor 24/7.

Create a fall monitoring system and protect patients by being alerted when a patient attempts to get out of bed, falls on a floor mat, or slides out of a seat. With the assistance of SPT-1, patients can be kept safe from harm, and family member can have a peace of mind that their loved ones will be protected.

## **Features**

- Compatible with bed, chair, and floor mat sensor systems.
- Detects when a patient gets up from the bed/chair pad, or steps on floor mat
- Connection type: standard handset jack (RJ-9)
- Dip switch for bed/chair pad or fall mat mode

- Wireless, small and lightweight
- Tamper protection against removal
- LED indicator
- Regular supervision signals to check system integrity
- Low battery detection

## **Specifications**

Frequency	433MHz / 869MHz
Power Source	3V, CR123 Lithium battery x1
Battery Life	3 years*
Operating Temperature	-10°C to 45°C (14°F to 113°F)
Operating Humidity	Up to 85% non-condensing
Dimensions	39mm x 78mm x 27mm

<sup>\*</sup> Note: Actual battery life may vary with product settings, operating environment, and usage.

